

# CCD Image Automatic Measuring system

## CCD Image Automatic Measuring system




### CCD Image Automatic Measuring system

#### Main Function:

CCD Image automatic measuring system for Vickers, Knoop and carburizing layer tests.

Equipped with CCD camera allow direct observation and measure of the indentation on the display. With its measuring software a single point measurement, random multipoint measurement, data statistical measurement, arbitrary two-point or multipoint measurement are possible and automatically display the depth, statistical calculations, conversions, display curves, data and images management, record and reports in Word or Excel documents.

Operating system: WIN7, WINXP



#### Standard Configuration:

1. Computer (Hard disk: 500G, memory: 2G, 19" LCD screen) 1 set
2. Ink jet printer 1 set
3. CCD camera 1 set
4. 1.5x adapter 1 pc
5. USB softdog 1 pc
6. Measuring software 1 pc



**DEVCO S.r.l.**

Via Marzabotto, 59 20037 Paderno Dugnano (MI)  
Tel +39 0283591153 - Fax. +39 0295441300  
[www.devcosrl.it](http://www.devcosrl.it) - e-mail: [info@devcosrl.it](mailto:info@devcosrl.it)

[www.echo-lab.it](http://www.echo-lab.it)