


## DIAMONDCUT-25

### Features

- Bench-top automatic diamond cutting machine, suitable for hard materials, especially suggested for artificial crystals
- Equipped with four different vices, allows the best cutting angle
- Easy to use with a compact body and safety protection
- High level of positioning accuracy
- Built-in cooling system with tank



Model		DIAMONDCUT-25
Cutting capacity		Max. ø 25mm
Cutting disc diameter		100 x 0,7 x 12,7mm
X axis travel		12mm, step 0.1mm
Noise level		<65dB
Working chamber LxW		196x104mm
Speed		50÷800rpm
Operating temperature		5÷30°C
Humidity non-condensing		≤85%
Motor power		600W
Power supply		230V ±10% - 50/60Hz - 1Ph
Tank capacity		0.5l
Machine	Dimensions	350x380x290mm
	Weight	25kg

### Standard accessories

Diamond cutting disc Ø100mm; Spanners