

SM-B20

Base for Stereo Microscope



Model	SM-B20
Base	With clip
	With black & white platen ø95mm
	Size 240x180x22mm
Drawtube	Height 280mm
	ø76mm
Focus distance	230mm
Power supply	110/220V 50/60Hz 1pH
Weight	4kg

SM-B21

Base for Stereo Microscope



Model	SM-B21
Base	With clip
	With glass platen ø95mm
	Size 240x180x22mm
Drawtube	Height 280mm
	ø76mm
Transmitted illumination	3W LED lamp
	Brightness control
Reflected illumination	3W LED lamp
	Brightness control
Focus distance	230mm
Power supply	110/220V 50/60Hz 1pH
Weight	6kg



DEVCO S.r.l.

Via Marzabotto, 59 20037 Paderno Dugnano (MI)
Tel +39 0283591153 - Fax. +39 0295441300
www.devcosrl.it - e-mail: info@devcosrl.it

www.echo-lab.it