

SM-B12**Base for Stereo Microscope**

Model	SM-B12
Base	With clip
	With glass platen ø95mm
	Size 285x238x22mm
Drawtube	Height 300mm
	ø76mm
Transmitted illumination	3W LED lamp
	Brightness control
Reflected illumination	3W LED lamp
	Brightness control
Focus distance	100mm
Power supply	110/220V 50/60Hz 1pH
Weight	3kg

SM-B14**Base for Stereo Microscope**

Model	SM-B14
Base	With clip
	With glass platen ø95mm
	Size 285x238x22mm
Pillar	Height 240mm
	ø32mm
Transmitted illumination	3W LED lamp
	Brightness control
Reflected illumination	3W LED lamp
	Brightness control
Arm type	SM-A1
Power supply	110/220V 50/60Hz 1pH
Weight	2kg



DEVCO S.r.l.

Via Marzabotto, 59 20037 Paderno Dugnano (MI)
Tel +39 0283591153 - Fax. +39 0295441300
www.devcosrl.it - e-mail: info@devcosrl.it

www.echo-lab.it