

VACUUM-Y**Features**

- Designed to be used with porous samples.
Allows the degassing of the material porosities and replacement with cold mounting resin
- Supplied with a fully enclosed transparent working chamber and diaphragm pump.



Model	VACUUM-Y
Vacuum chamber dimensions	ø 180x314mm
Vacuum level	0.087MPa
Noise level	66dB
Operating temperature	7÷40°C
Humidity non-condensing	65%
Motor power	180W
Power supply	230V ±10% - 50/60Hz - 1Ph
Dimensions	480x400x400mm
Weight	40kg

Standard accessories

Set of plastic pipes; epoxy resin; hardener; diaphragm pump



DEVCO S.r.l.

Via Marzabotto, 59 20037 Paderno Dugnano (MI)
Tel +39 0283591153 - Fax. +39 0295441300
www.devcosrl.it - e-mail: info@devcosrl.it

www.echo-lab.it