

PNM 100

Manual notching machine for Charpy and Izod plastic tests specimen

FOR PLASTIC SAMPLE



Cutting tool:

- > Type A in conformity to ISO 179, ISO 180, ASTM D256, ASTM D6110
- > Type B in conformity to ISO 179, ISO 180

Included accessories

Knife to choose between Type A or Type B
Fixture set

Accessories

| Model | Description |
|-------|--|
| BR01A | Knife type A compliance to ISO 179, ISO 180, ASTM D256, ASTM D6110 |
| BR01B | Knife type B compliance to ISO 179, ISO 180 |

| Model | PNM 100 | |
|---|-----------|--------------|
| Motor speed | rpm | 1440 |
| Feeding stroke (specimen max thickness) | mm | 25 |
| Power supply | V-Hz-W-ph | 230-50-180-1 |
| Dimensions | mm | 326x376x460 |
| Weight | kg | 30 |