

ADDITIONAL INSTRUMENTS



Model	Type	Stroke/min	Power	Power supply	Standard accessories	Weight	Recommended for:
EG 710	Engraving pen	7100	0,1 kW	230V - 50Hz	2 set of tips	0,3kg	Samples marking

Model	Type	Heating power	Power supply	Recommended for:
DD-1	Drying device	2000W, 900W, ordinary temperature	230V - 50Hz	Fast sample drying

Model	Description	Recommended for:
SHM-100	Sample holder for round specimens (diameter on request). Mechanical type, with adjustable samples height	Easy handling of specimens

Model	Description	Recommended for:
D-1	Automatic dispenser for diamond/alumina/silica products	Automatic distribution of products during grinding/polishing processes