

ZETAPRESS M2D

Features

- > Automatic mounting press with large size LCD touch screen, microprocessor controlled for setting all process parameters such as, automatic pre-heating, heating, working pressure, cooling time
- > One cylinder with intermediate piston which makes possible to carry out 2 mounting operations at same time
- > Different mold sizes available (\varnothing 25mm, \varnothing 32mm, \varnothing 40mm, \varnothing 50mm)
- > Automatic water cooling system



Model	ZETAPRESS M2D
Mold diameter	\varnothing 25, \varnothing 30, \varnothing 40, \varnothing 50mm
Mounting quantity per time	2
Moulding temperature	50÷300
Mounting time	<15min
Heating time to 100°C	8min
Heating time to 140°C	10min
Heating time to 180°C	12min
Heating power	1000W
Cooling time	3÷5min
Pressure	0÷72MPa
Noise level	65dB
Operating temperature	0÷30°C
Humidity non-condensing	<85%
Water inlet	8mm
Water outlet	13mm
Power supply	230V \pm 10% - 50/60Hz - 1Ph
Dimensions	630x520x520mm
Weight	120

Standard accessories

To be selected from: \varnothing 25mm, \varnothing 30mm, \varnothing 40mm, \varnothing 50mm; mounting powder; filler; inlet water pipe; outlet water pipe; hose clamp