

Microscopes
accessories

M-TCS

LS-CLS

Temperature controlled stage

Cold Light Source



Specifications	M-TCS
Pillar height	240mm
Diameter	32mm
Base size	380x260x20mm
Power supply	230V ±10% - 50/60Hz -1Ph - 50W
Temperature control range	20°C (room Temperature) to 50°C
Temperature setting Div	0.1°C
Variance between setting and displaying Temp	<+/-1°C
Variance between actual stage temperature and displaying Temp	<+/-1°C
Stage dimensions	150x130x5.8mm

Specifications	LS-CLS
Power Supplier	Power Switch 110V/ 220V - T-H Lamp with Reflector 20V/ 150W - Adjustable Brightness - Multi-coated IR Cutoff Filter for Heat Insulation - Fan Cooling
Dimensions	260x120x150mm
Weight	3.3 kg

LS-C40WL

Cold 40W Light Source



LS-2L3Wpf

Post Fixing Light Source, 3W LED 2



Specifications	LS-C40WL
Light Source	Power LED 40W
Brightness	Continuously adjustable
Working life	More than 20000 hours
Power supply	230V ±10% - 50/60Hz - 1Ph - 50W
Light Guide	Single, Rigid Flexional Tube, Length 500mm - Dual, Rigid Flexional Tube, Length 500mm - Ring, Flexible Tube, Length 1000mm, Size for connecting Microscope Ø60mm
Kohler	Kohler Illumination Attachment
Power	Wide range voltage 100-240V

Specifications	LS-2L3Wpf
Luminous Flux	240LM
LED Power	3W
LED Quantity	2
Fiber length	450mm
Light Intensity Control	Continuously
Working Temperature	≤50°C
Input Voltage	90-265V
Size	140mmx98mmx52mm

LS-C5WL

Cleaning kit

Cold 5W LED Light Source 5600K

Complete kit for cleaning of eyepieces, objectives and microscopes components.



Specifications	LS-C5WL
Light Source	Power LED 5W
Brightness	Continuously adjustable
Power	85V- 240V,50- 60Hz
Optical Interface	With a 5/8" Fiber Optic Cable adapter with metal housing
Color Temperature	5600K
Dimensions	200x85x110mm



DEVCO S.r.l.

Via Marzabotto, 59 20037 Paderno Dugnano (MI)
Tel +39 0283591153 - Fax. +39 0295441300
www.devcosrl.it - e-mail: info@devcosrl.it

www.echo-lab.it